

AORTIC ENDOGRAFTS

AAA MAIN BODY GRAFTS

Company	Device	Main Body		Iliac Leg(s)	
		Length (cm)	Diameter (mm)	Length (cm)	Diameter (mm)
Cook Medical	Zenith Flex	11.2, 12.5,* 12.6, 14.1, 14.3,* 15.5, 16.1,* 17, 17.9*	22, 24, 26, 28, 30, 32, 36	3.7,† 3.9, 5.4,‡ 5.6, 7.1,‡ 7.3, 8.8,‡ 9, 10.5,‡ 10.7,§ 12.2,‡ 12.4§	8, 10, 12, 14, 16, 18, 20, 22, 24
Endologix	Powerlink	12, 13.5, 14, 15.5, unibody main sizes 8, 10	25, 28	4, 5.5	16
Medtronic Vascular	AneuRx AAAAdvantage	13.5, 16.5	20, 22, 24, 26, 28	8.5, 11.5, 13.5	12, 13, 14, 15, 16, 18, 20 (20, 22, 24 distal flare)
	Talent Abdominal	14, 15.5, 17	22, 24, 26, 28, 30, 32, 34, 36	7.5, 9, 10.5	8, 10, 12, 14, 16, 18, 20, 22, 24
W. L. Gore & Associates	Gore Excluder AAA Endoprosthesis	12, 13,¶ 14, 15,¶ 16, 17,¶ 18	23, 26, 28.5, 31†	9.5,‡ 10, 11.5,‡ 12, 13.5,‡ 14	12, 14.5, 16, 18, 20

*Lengths for 36-mm diameter grafts only.

†Sizes available outside US.

‡Lengths for 8- and 10-mm diameter grafts only.

§Lengths for 12- and 14-mm diameter grafts only.

¶Lengths for 31-mm diameter grafts only.

‡Lengths for 16-, 18-, and 20-mm diameter grafts only.

AAA ANCILLARY COMPONENTS

Company	Product	Main Body Extensions		Iliac Leg Extensions
		Length (cm)	Diameter (mm)	Length (cm)
Cook Medical	Zenith	3.9, 5,* 5.8, 7.3*	22, 24, 26, 28, 30, 32, 36	5.5
	Renu	4.3, 5.4,* 6.2, 7.7,* 8.1, 10*		N/A
Endologix	Powerlink	5.5, 7.5, 8, 9.5, 10	25, 28, 34	5.5, 6.5, 8.8
Medtronic Vascular	AneuRx AAAAdvantage	4	20, 22, 24, 26, 28	5.5, 8.5
	Talent Abdominal	2.6, 2.8, 2.9, 3.0	22, 24, 26, 28, 30, 32, 34, 36	7.4, 7.5, 7.9, 8, 8.1, 14
W. L. Gore & Associates	Gore Excluder AAA Endoprosthesis	3.3, 4.5†	23, 26, 28.5, 32†	7, 9.5,‡ 11.5,‡ 13.5‡

*Lengths for 36-mm diameter grafts only.

†Sizes available outside US.

‡Lengths for 8- and 10-mm diameter grafts only.

§Lengths for 12- and 14-mm diameter grafts only.

¶Lengths for 31-mm diameter grafts only.

‡Lengths for 16-, 18-, and 20-mm diameter grafts only.

AORTIC ENDOGRAFTS

	Main Body Delivery System Profile			Fixation Location	Stent Expansion	Stent Material	Graft Material	Markets Available
	Device Outer Diameter (OD) (F)	Sheath Required for Delivery?	Delivery Sheath OD (F)					
	20, 23, 26	No	N/A	Suprarenal	Self-expanding	Stainless steel	Woven polyester	US, Europe, Asia, Australia, Canada, South America, Middle East, Africa, Japan
	21, 22 [†]	No	N/A	Infrarenal	Self-expanding	Cobalt-chromium alloy	High-density ePTFE	US, Europe, Latin America, Africa, Japan
	21	No	N/A	Infrarenal	Self-expanding	Nitinol	Woven polyester	US, Latin America, Asia
	22, 24			Suprarenal				US, Europe, Middle East, Africa, Latin America, Asia, Australia, Canada, India
	18, 20 [†]	Yes	20.4, 22.8 [†]	Infrarenal	Self-expanding	Nitinol	ePTFE	US, Canada, Europe, Africa, Asia, Australia, Latin America, Japan, Middle East, South America, India

	Aorto-uni-iliac/Converter		Occluders
Diameter (mm)	Length (cm)	Proximal Diameter (mm)	Diameter (mm)
8, 10, 12, 14, 16, 18, 20, 22, 24	8, 8.2*	24, 28, 32, 36	14, 16, 20, 24
N/A	9.1, 9.9,* 10.8, 11.3, 11.6,* 12.5, 12.7,* 13.0, 13.3,* 14.4,* 14.7, 16.1*	22, 24, 26, 28, 30, 32, 36	N/A
16, 20, 25	N/A	N/A	N/A
12, 13, 14, 15, 16, 18, 20 (20, 22, 24 distal flare)	N/A	N/A	N/A
8, 10, 12, 14, 16, 18, 20, 22, 24			
10, 12, 14.5, 16, 18, 20	N/A	N/A	N/A