

ENDOVASCULAR TREATMENT AND MEASUREMENT DEVICE

Company Name	Product Name	Measurement Range (mm)	Measurement Balloon Length (mm)	PTCA Balloon Diameter (mm)	PTCA Balloon Length (mm)	PTCA Nominal Pressure (atm)	Rated Burst Pressure (atm)	Catheter Usable Length (cm)	Maximum Guidewire (inch)	Minimum Guide Catheter ID (inch)
Neovasc Inc.	Metricath Gemini	2-4	6	2.5, 3, 3.5, 4	15	6	16	140	.014 (.356 mm)	.064 (6 F)

Comments

Treatment balloon combined with a low-pressure measurement balloon that quickly and easily determines precise vessel and in-stent diameter and cross-sectional area

FDA Indicated Use

510(k) cleared for use in the renal, iliac, and femoral arteries