

OTHER DEVICES

ATHERECTOMY DEVICES

Company Name	Product Name	Minimum Vessel Diameter (mm)	Sheath Compatibility (F)	Crossing Profile (inch)	Working Length (cm)	Tip Length (cm)	
Cardiovascular Systems, Inc. (CSI)	Diamondback 360° 1.25-mm Classic Crown	1.5	6	.069 (1.75 mm)	135	1	
	Diamondback 360° 1.5-mm Classic Crown					1.5	
	Diamondback 360° 1.75-mm Classic Crown						
	Diamondback 360° 2-mm Classic Crown					.079 (2 mm)	2
	Diamondback 360° 2.25-mm Classic Crown			7		.089 (2.25 mm)	
	Diamondback 360° 1.5-mm Solid Crown			6		.069 (1.75 mm)	1
	Diamondback 360° 1.75-mm Solid Crown						1.5
	Diamondback 360° 2-mm Solid Crown					.079 (2 mm)	
	Diamondback 360° 2.25-mm Solid Crown			7		.089 (2.25 mm)	2
ev3 Inc.	SilverHawk LS-C	4.5–7	8	.105 (2.7 mm)	107	6	
	SilverHawk LS-M		7, 8				
	SilverHawk LS-F		8				
	SilverHawk LX-M	4.5–6.5	7, 8		113	9	
	SilverHawk MS-M	3.5–5			110	6	
	SilverHawk MS-F		8		107		
	SilverHawk SX	3–3.5	7		.095 (2.4 mm)	136	4.3
	SilverHawk SXL					136	7.2
	SilverHawk SS+				.09 (2.3 mm)	135	2.6
	SilverHawk EXL	2–3	6		.08 (2 mm)		6
	SilverHawk ES+	2–2.5			.075 (1.9 mm)		2.2
	SilverHawk DS	1.5–2			.073 (1.9 mm); .062 (1.57 mm) without wire		2.6
	Pathway Medical Technologies	Jetstream	3		8	.083 (2.1 mm)	135