

# OTHER DEVICES

TARGETED RENAL THERAPY						
Company Name	Product Name	Sheath Size (F)	System Lengths (cm)	Maximum Guiding Catheter Size (F)	Bifurcated Infusion Branch Profile (inch)	Indicated Use
FlowMedica, Inc.	Benephit PV Infusion System	5	40, 105, 140	N/A	.035	Physician-specified agents to the peripheral vasculature including but not limited to the renal arteries
	Benephit XT Infusion System		35, 105			
	Benephit CV Infusion System	8	30, 35, 40, 45	6		