

PTA BALLOONS

Company Name	Balloon Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Size (F)	Sheath Sizing Size (F)	Shaft Length (cm)	Balloon Diameters (mm)	Balloon Lengths (cm)
Abbott Vascular	Fox Plus	8	up to 18	.038	5-7	5-7	80, 135	3-12	20-80
	Fox SV	6	up to 22	.018	4	4	90, 135	2, 2.5, 3, 4, 5	15, 20, 30, 40, 60
AngioDynamics	WorkHorse I/II	10	14, 16, 18	.035	5	5-8 mm through a 6 F; 9-10 mm through a 7 F	50, 75, 100	5, 6, 7, 8, 9, 10	2, 4, 6
	Profiler	7	12, 14, 16	.035	5, 5.7	4-8 mm X 4 mm through a 5 F; 8 mm X 6-10 mm through a 6 F	40, 75, 120	4, 5, 6, 7, 8, 9, 10	2, 4, 6, 8
Bard Peripheral Vascular, Inc.	Conquest	8	up to 30	.035	Not provided	6-8	50, 75	5, 6, 7, 8, 9, 10, 12	2, 3, 4, 6, 8
	Atlas	4, 6	up to 18	.035	Not provided	7-12	75, 120	12, 14, 16, 18, 20, 22, 24, 26	2, 4, 6
	Centurion	10	up to 20	.035	Not provided	6-8	50, 75, 120	4, 5, 6, 7, 8, 9, 10	2, 4
	Optiplast XT	5, 6, 8	up to 15	.035	Not provided	5-8	50, 75, 100, 120	3, 4, 5, 6, 7, 8, 9, 10, 12	1.5, 2, 3, 4, 6, 8, 10
	Ultraverse	8	up to 15	.018	Not provided	4, 5	100, 120	2, 3, 4, 5, 6, 7	2, 4, 6, 10
Boston Scientific Corporation	Symmetry	6	15	.018	4, 5	4 (1.5-4 X 4 cm), 5 (4 X 10 cm-6 X 4 cm)	90, 135	1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6	2, 4, 10
	Symmetry Stiff Shaft	6	15	.018	4, 5	5 (1.5-4 X 4 cm), 5 (4 X 10 cm-6 X 4 cm)	90, 135, 150	1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 6	2, 4, 10
	Talon	7	15 (4-6 mm), 12 (7 mm)	.018	5, 6	5 (4-6 mm), 6 (7 mm)	90, 135	4, 5, 6, 7	1.5, 2, 4
	Gazelle	7 (2-4 mm), 8 (4.5-6 mm)	14	.018	4	4	90, 135	2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6	2
	Ultra-Soft SV Monorail	6	14	.018	4, 5	4 (4-5 mm), 5 (5.5-7 mm)	90, 150	4, 4.5, 5, 5.5, 6, 6.5, 7	1.5, 2
	XXL	Not provided	8 (12, 14 mm), 5 (16, 18 mm)	.035	7, 8	8 (12, 14 mm), 8 (16, 18 mm)	75, 120 mm	12, 14, 16, 18	2, 4, 6

PTA BALLOONS

Company Name	Balloon Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Size (F)	Sheath Sizing Size (F)	Shaft Length (cm)	Balloon Diameters (mm)	Balloon Lengths (cm)
Boston Scientific Corporation	Blue Max 20	6	20 (4-8 X 4 cm), 17 (8 X 8-10 cm)	.035	6, 7, 8	6 (3-6 mm), 7 (7-8 mm), 8 (9-10 mm)	40, 75, 120	4-10	2, 3, 4, 8, 10
	Ultra-thin Diamond	6	15 (3-8 mm), 12 (9-10 mm), 10 (12 mm)	.035	5, 6, 7	15 (3-6 mm), 12 (7-8 mm), 10 (9-12 mm)	40, 75, 120, 135, 150	3-12	1.5, 2, 3, 4, 6, 10
	Ultra-thin SDS	10	12	.035	5, 6, 7	5 (up to 6 X 4), 6 (up to 8 X 4), 7 (8 X 6-10 cm)	50, 75, 90, 135	4-10	1.5, 2, 3, 4, 6, 8
	Synergy	6 (3-10 mm), 4 (12 mm)	12 (3-8 mm), 10 (9-10 mm), 6 (12 mm)	.035	5, 6, 7	5 (up to 6 X 4), 6 (up to 8 X 4), 7 (8 X 6-12 cm)	50, 75, 90, 135, 150	3-12	1.5, 2, 3, 4, 6, 8, 10
Cook Incorporated	ATB Advance	5	15	.035	5	Not provided	40, 80, 120	4, 5, 6, 7, 8, 9, 10	2, 3, 4, 5, 6, 7, 8
Cordis Endovascular, a division of Cordis Corporation	Opta Pro PTA Dilatation Catheter	6, 4	10, 6	.035	5	5-7	80, 110, 135	3-10, 12	1, 1.5, 2-4, 6, 8, 10
	Powerflex P3 PTA Dilatation Catheter	10, 6	15, 8	.035	5	5-8	40, 65, 80, 110, 135	4-10, 12	1, 1.5, 2-4, 6, 8, 10
	Powerflex Extreme PTA Dilatation Catheter	15, 12	20, 17	.035	5	5-9	40, 80, 120	4-10	2, 4, 6
	Maxi LD PTA Dilatation Catheter	4, 3.5	6, 5	.035	7	8-11	80, 110	Iliac: 14, 15 Esophageal: 16, 18, 20, 22, 25	2, 4, 6, 8
	Slalom PTA Dilatation Catheter	10, 8, 6	14, 12, 10	.018	3.7, 4	4-6	80, 120, 135	3-8	2, 4
	Savvy PTA Dilatation Catheter	8, 6	10	.018	3-3.5	4-5	80, 120, 150	2-6 in .5-mm increments	2-4, 6 and 10

PTA BALLOONS

Company Name	Balloon Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Size (F)	Sheath Sizing Size (F)	Shaft Length (cm)	Balloon Diameters (mm)	Balloon Lengths (cm)	
Cordis Endovascular, a division of Cordis Corporation	Aviator PTA Dilatation Catheter	10, 8	14, 12	.014	33/33-29, 33/33, 36/33, 33/35, 33/35, 36/35	4-6	135	4-7 in .5-mm increments	1.5, 2, 3, 4	
ev3 / Invatec	Amphirion Deep OTW PTA Balloon Catheter	7	14	.014	2.9	4	120, 150	1.5, 2, 2.5, 3, 3.5, 4	2, 4, 8, 12	
	Sailor Plus OTW PTA Balloon Catheter	7	up to 14	.035	5	5-7	40, 80, 130	3, 4, 5, 6, 7, 8, 9, 10, 12	2, 3, 4, 6, 8, 12	
	Submarine Plus OTW PTA Balloon Catheter	7	up to 14	.018	3.7	4-5	90, 130	2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6, 7	2, 4, 6, 8, 12	
Guidant Corporation	Viatrac 14 Plus Peripheral Dilatation Catheter	8	14	.014	6	4	80, 135	4	1.5, 2, 3, 4	
						5		4.5		
								5, 5.5, 6, 6.5, 7		
	Agiltrac 035 Peripheral Dilatation Catheter	8	14	.035	7	5	55, 80, 135	4	2, 4, 6, 8, 10	
								5	2, 4, 6, 8, 10	
								6	2, 4, 6, 8	
						8		6	10	
						7		5	2	
						8		6	4, 6, 8, 10	
								8	2, 4, 6, 8	
								9	2, 4, 6	
								10	2, 4	
								9	6	
						6		8	2, 4	
								14	2, 4	
Agiltrac 018 Peripheral Dilatation Catheter	8	14	.018	6	5	80, 135	4, 5, 6, 7	2		
					7		8			
					8		9, 10			
							6	5	4, 5, 6	3, 4
					7		6	7, 8		
					8			9, 10		
							8	9	12	6
							14	7	7	
							8	6	8	
							12	9	9	
							12	7	10	
					Medtronic Vascular		Sprinter	8, 6 (1.50 mm)	14, 12 (1.50 mm)	n/a
NC Stormer	10	18	n/a	5		n/a	138	2, 2.25, 2.50, 2.75, 3, 3.25, 3.50, 3.75, 4	6, 11, 16, 21, 26, 31	

EMBOLECTOMY BALLOON/CATHETER						
Company Name	Balloon Name	Size (F)	Sheath Sizing Size (F)	Shaft Length (cm)	Balloon Diameters (mm)	Balloon Lengths (cm)
Atrium Medical Corporation	ClearWay TBI Catheter	6.5	8	40, 80	6, 7, 8	1
Edwards Lifesciences	Edwards Fogarty Adherent Clot Catheter	4	6 or 7 (dependent upon sheath)	80	6	n/a
		Not provided				
	Edwards Fogarty Embolectomy Catheter	2	3.5	60	4	Not provided
		3	4	40, 80	5	
		4	5	40, 80	9	
		5	6	80	11	
		6	7	80	13	
		7	8	80	14	
		Edwards Over-the-Wire Fogarty Embolectomy Catheter	3	4.5	40, 80	
	4		5	40, 80	9	
	5.5		6	40, 80	11	
	6		7	40, 80	13	
	7		8	80	14	
	LeMaitre Vascular, Inc.	LeMaitre Single Lumen Embolectomy Catheter	2	Not provided	40, 60	4.5
3			40, 80		8	1.4
4			40, 80		10.5	1.7
5			40, 80		13	2
6			40, 80		13.5	2.25
7			40, 80		14	2.4
NovaSil Latex Free Single Lumen Embolectomy Catheter			2		Not provided	40, 60
		3	40, 80	6		1.4
		4	40, 80	9		1.7
		5	40, 80	11		2
		6	40, 80	13		2.25
LeMaitre Over-the-Wire Embolectomy Catheter		3	4.5	40, 80	6	1.4
		4	5	40, 80	10	1.7
		5	6	40, 80	12	2
		5+	6	40, 80	12	2
		6	7	40, 80	13	2.25
		7	8	40, 80	14	2.4

PTA BALLOONS

SPECIALTY BALLOONS									
Company Name	Balloon Name	Nominal (atm)	Rated Burst (atm)	Endhole (in)	French Size	Sheath Sizing French Size	Shaft Length (cm)	Balloon Diameters (mm)	Balloon Lengths (cm)
Boston Scientific Corporation	PolarCath Cryoplasty Dilatation System	8	n/a	.035	6, 7	6 (4-5 mm), 7 (6-7 mm), 8 (8 mm)	80, 120	4-8	2, 4, 6
				.014	5	5	135	2.5, 3, 3.5, 4	2, 4, 6
	Peripheral Cutting Balloon	6	10	.018	6, 7	6 (5, 6 mm), 7 (7, 8 mm)	50, 90, 135	5-8	2
				.014	6	6	140	2, 2.5, 3, 3.5, 4	1.5