

PTA BALLOONS

Company Name	Balloon Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Size (F)	Sheath Sizing Size (F)	Shaft Length (cm)	Balloon Diameter (mm)	Balloon Length (cm)																	
Abbott Vascular	Viatrac 14 Plus Peripheral Dilatation Catheter	8	14	.014	6	4	80, 135	4, 4.5	1.5, 2, 3, 4																	
						5		5, 5.5, 6, 6.5, 7																		
	Agiltrac .035 Peripheral Dilatation Catheter	8	14	.035	7	5	55, 80, 135	4	2, 4, 6, 8, 10																	
						6		2, 4, 6, 8, 10																		
						8		2, 4, 6, 8																		
						10		10																		
						12		2																		
						14		4, 6, 8, 10																		
						16		2, 4, 6, 8																		
						18		2, 4, 6																		
						20		2, 4																		
						22		6																		
						24		2, 4																		
						26		2, 4																		
	Agiltrac .018 Peripheral Dilatation Catheter	8	14	.018	6	5	80, 135	4, 5, 6, 7	2																	
						7		8																		
						8		9, 10																		
						9		4, 5, 6		3, 4																
						10		7, 8																		
						11		9, 10																		
12						12																				
14						8		6																		
16						4, 5, 6, 7																				
18						8																				
20						9																				
Fox Plus						8		Up to 18	.035	5-7	5-7	80, 135	3-12	20-80												
	Fox sv	6	Up to 22	.018	4		4								90, 135	2, 2.5, 3, 4, 5	15, 20, 30, 40, 60, 80, 100, 120									
																		WorkHorse I/II	10	14, 16, 18	.035	5	5-8 mm through a 6-F, 9-10 mm through a 7-F	50, 75, 100	5, 6, 7, 8, 9, 10	2, 4, 6

PTA BALLOONS

Company Name	Balloon Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Size (F)	Sheath Sizing Size (F)	Shaft Length (cm)	Balloon Diameter (mm)	Balloon Length (cm)
Bard Peripheral Vascular, Inc.	Conquest	8	Up to 30	.035	Not provided	6-8	50, 75	5, 6, 7, 8, 9, 10, 12	2, 3, 4, 6, 8
	Atlas	4, 6	Up to 18	.035	Not provided	7-12	75, 120	12, 14, 16, 18, 20, 22, 24, 26	2, 4, 6
	Centurion	10	Up to 20	.035	Not provided	6-8	50, 75, 120	4, 5, 6, 7, 8, 9, 10	2, 4
	Optiplast XT	5, 6, 8	Up to 15	.035	Not provided	5-8	50, 75, 100, 120	3, 4, 5, 6, 7, 8, 9, 10, 12	1.5, 2, 3, 4, 6, 8, 10
	Ultraverse	8	Up to 15	.018	Not provided	4, 5	100, 120	2, 3, 4, 5, 6, 7	2, 4, 6, 10
Biotronik Vascular Inc.	Pheron	7	15	.035	5, 6	5-6	80, 130	3, 4, 5, 6, 7, 8, 9, 10	2, 4, 6, 8, 10
B. Braun	Mini Ghost	up to 10	Up to 20	.018	3.5	4	40, 80, 120, 150	2-6	2, 3, 4, 6, 8, 10
	Z-Med	up to 6	Up to 15	.035	5-11	5-16	80, 100, 120	2-40	2, 3, 4, 5, 6
	Impact	up to 6	Up to 10	.035	6-9	8-14	75, 120	12-25	2, 4
	CoEfficient	up to 6	Up to 13	.018	3.5	4-7	75, 135	4-12	2, 3, 4
Boston Scientific Corporation	Symmetry Balloon Catheter with Glidex Hydrophilic Coating	6	15	.018	4, 5	4 (1.5 mm-4 mm X 4 cm), 5 (4 mm X 10 cm-6 mm X 4 cm)	90, 135	1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6	2, 4, 10
	Symmetry Stiff Shaft with Glidex Hydrophilic Coating Balloon Catheter	6	15	.018	4, 5	4 (1.5 mm-4 mm X 4 cm), 5 (4 mm x 10 cm-6 mm X 4 cm)	90, 135, 150	1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 6	2, 4, 10
	Talon	7	15 (4-6 mm), 12 (7 mm)	.018	5, 6	5 (4-6 mm), 6 (7 mm)	90, 135	4, 5, 6, 7	1.5, 2, 4
	Gazelle	7 (2-4 mm), 8 (4.5-6 mm)	14	.018	4	4	90, 135	2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6	2
	Ultra-Soft SV Monorail	6	14	.018	4, 5	4 (4-5 mm), 5 (5.5-7 mm)	90, 150	4, 4.5, 5, 5.5, 6, 6.5, 7	1.5, 2

PTA BALLOONS

Company Name	Balloon Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Size (F)	Sheath Sizing Size (F)	Shaft Length (cm)	Balloon Diameter (mm)	Balloon Length (cm)
Boston Scientific Corporation	XXL Balloon Catheter	Not provided	8 (12, 14 mm), 5 (16, 18 mm)	.035	7, 8	7 (12, 14 mm), 8 (16, 18 mm)	75, 120	12, 14, 16, 18	2, 4, 6
	Blue Max 20 Balloon Catheter	6	20 (4-8 mm X 4 cm), 17 (8 mm X 8 cm-10 mm)	.035	6, 7, 8	6 (4-6 mm), 7 (7-8 mm), 8 (9-10 mm)	40, 75, 120	4, 5, 6, 7, 8, 9, 10	2, 3, 4, 8, 10
	Ultra-thin Diamond Balloon Catheter	6	15 (3-8 mm), 12 (9-10 mm), 10 (12 mm)	.035	5, 6, 7	5 (3-6 mm), 6 (7-8 mm), 7 (9-12 mm)	40, 75, 120, 135, 150	3, 4, 5, 6, 7, 8, 9, 10, 12	1.5, 2, 3, 4, 6, 8, 10
	Ultra-thin SDS Balloon Catheter	6-8	12	.035	5, 6, 7	5 (up to 6 X 4), 6 (up to 8 X 4), 7 (9-10 mm)	50, 75, 90, 135	4, 5, 6, 7, 8, 9, 10	1.5, 2, 3, 4, 6, 8
	Synergy Balloon Catheter	6 (3-10 mm), 4 (12 mm)	12 (3-8 mm), 10 (9-10 mm), 6 (12 mm)	.035	5, 6, 7	5 (up to 6 X 4), 6 (up to 8 X 4), 7 (9-12 mm)	50, 75, 90, 135, 150	3, 4, 5, 6, 7, 8, 9, 10, 12	1.5, 2, 3, 4, 6, 8, 10
	Sterling Monorail PTA Balloon Dilatation Catheter	6	14 (3-7 mm), 12 (8 mm)	.018, .014	4, 5	4 (up to 8 X 4), 5 (up to 8 X 6)	80, 135	3, 3.5, 4, 4.5, 5, 5.5, 6, 6.5, 7, 8	1, 1.5, 2, 3, 4, 6
	Sterling Over-the-Wire Balloon Dilatation Catheter	6	14 (4-7 mm), 12 (8 mm), 10 (9-10 mm)	.018, .014	4, 5, 6	4 (up to 8 X 4), 5 (up to 10 X 4), 6 (up to 10 X 8)	40, 80, 135	4, 5, 6, 7, 8, 9, 10	2, 3, 4, 6, 8, 10
Cook Incorporated	ATB Advance	5	15	.035	5	5-7	40, 80, 120	4, 5, 6, 7, 8, 9, 10, 12, 14	2, 3, 4, 5, 6, 7, 8
Cordis Endovascular Division of Cordis Corporation	Opta Pro PTA Dilatation Catheter	6, 4 (12 mm)	10, 6 (12 mm)	.035	5	5-7	80, 110, 135	3-10, 12	1, 1.5, 2-4, 6, 8, 10
	Powerflex P3 PTA Dilatation Catheter	10, 6 (12 mm)	15, 8 (12 mm)	.035	5	5-8	40, 65, 80, 110, 135	4-10, 12	1, 1.5, 2-4, 6, 8, 10

PTA BALLOONS

Company Name	Balloon Name	Nominal (atm)	Rated Burst (atm)	Endhole (inch)	Size (F)	Sheath Sizing Size (F)	Shaft Length (cm)	Balloon Diameter (mm)	Balloon Length (cm)
Cordis Endovascular Division of Cordis Corporation	Powerflex Extreme PTA Dilatation Catheter	15, 12 (10 mm)	20, 17 (9-10 mm)	.035	5	5-9	40, 80, 120	4-10	2, 4, 6
	Maxi LD PTA Dilatation Catheter	4, 3.5 (22, 25 mm)	6, 5 (22, 25 mm)	.035	7	8-11	80, 110	Iliac: 14, 15 Esophageal: 16, 18, 20, 22, 25	2, 4, 6, 8
	Slalom PTA Dilatation Catheter	10 (3-6 mm) 8 (6-7 mm) 6 (8 mm)	14 (3-5 mm) 12 (6-7 mm) 10 (8 mm)	.018	3.7, 4	4-6	80, 120, 135	3-8	2, 4
	Savvy PTA Dilatation Catheter	8(24mm) 6(46mm)	10	.018	3.5-3	4-5	80, 120, 150	2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6	2-4, 6, 10
	Aviator PTA Dilatation Catheter	10 (4-6.5 mm), 8 (7 mm)	14 (4-6.5 mm), 12 (7 mm)	.014	3.3/3.3-2.9, 3.3/3.3, 3.6/3.3, 3.3/3.5, 3.3/3.5, 3.6/3.5	4-6	135	4, 4.5, 5, 5.5, 6, 6.5, 7	1.5, 2, 3, 4
ev3 / Invatec	Amphirion Deep OTW PTA Balloon Catheter	7	14	.014	2.9	4	120, 150	1.5, 2, 2.5, 3, 3.5, 4	2, 4, 8, 12
	Admiral Xtreme OTW PTA Balloon Catheter	8, 6 (10, 12 mm)	Up to 14	.035	5	5 (3-7 mm), 6 (8-10 mm), 7 (12 mm)	80, 130	3, 4, 5, 6, 7, 8, 9, 10, 12	2, 4, 6, 8, 12
	Sailor Plus OTW PTA Balloon Catheter	7	Up to 14	.035	5	5 (3-7 mm), 6 (8 mm), 7 (9-12 mm)	40, 80, 130	3, 4, 5, 6, 7, 8, 9, 10, 12	2, 3, 4, 6, 8, 12
	Submarine Plus OTW PTA Balloon Catheter	7	Up to 14	.018	3.7	4 (2-5 X 40, 5.5), 5 (5 X 60-7)	90, 130	2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6, 7	2, 4, 6, 8, 12
Medtronic Vascular	Sprinter	8, 6 (1.5 mm)	14, 12 (1.5 mm)	n/a	5	n/a	138	1.5, 2, 2.25, 2.5, 2.75, 3, 3.25, 3.50, 3.75, 4	6, 10, 12, 15, 20, 25, 30
	NC Stormer	10	18	n/a	5	n/a	138	2, 2.25, 2.5, 2.75, 3, 3.25, 3.5, 3.75, 4	6, 11, 16, 21, 26, 31