

EMBOLECTOMY BALLOON/CATHETER

Company Name	Balloon Name	Size (F)	Sheath Sizing	Shaft Length (cm)	Balloon Diameter (mm)	Balloon Length (cm)	
Arrow International	Single Lumen	3	5	40	5	Not provided.	
		3	5	80	5		
		4	6	40	9		
		5	7	80	11		
		5	7	40	11		
	Double Lumen OTW	4	6	40	9		
		4	6	80	9		
		5	7	40	11		
5		7	80	11			
Atrium Medical Corporation	ClearWay TBI Catheter	6.5	8	40, 80	6, 7, 8	1	
Edwards Lifesciences	Edwards Fogarty Adherent Clot Catheter	4	6 or 7 (dependent upon sheath)	80	6	n/a	
	Edwards Fogarty Embolectomy Catheter	2	3.5	60	4	Not provided.	
		3	4	40, 80	5		
		4	5	40, 80, 100	9		
		5	6	80	11		
		6	7		13		
		7	8		14		
	Edwards Over-the-Wire Fogarty Embolectomy Catheter	3	4.5	40, 80, 120	5	Not provided.	
		4	5	40, 80	9		
		5.5	6		11		
		6	7		13		
		7	8	80	14		
	LeMaitre Vascular, Inc.	LeMaitre Single Lumen Embolectomy Catheter	2	Not provided.	40, 60	4.5	1.1
			3		40, 80	8	1.4
4					10.5	1.7	
5					13	2	
6					13.5	2.25	
7					14	2.4	
NovaSil Latex Free Single Lumen Embolectomy Catheter			2		Not provided.	40, 60	4
		3	40, 80	6		1.4	
		4		9		1.7	
		5		11		2	
		6		13		2.25	
		7		14		2.4	
LeMaitre Over-the-Wire Embolectomy Catheter		3	4.5	40, 80	6	1.4	
		4	5		10	1.7	
		5	6		12	2	
		6	7		13	2.25	
		7	8		14	2.4	