

OTHER DEVICES

CTO DEVICES

Company Name	Product Name	Catheter Size (F)	Wire Size (inch)	Working Length (cm)	Minimum Guiding Catheter ID (inch)	Comments
Cordis Endovascular Division of Cordis Corporation	Frontrunner XP CTO Catheter	3.1	n/a	90, 140	.067	2.3 mm, actuating distal tip.
	Frontrunner Micro Guide Catheter	4.5	n/a	82, 132	.078	Support catheter for use with 90-cm or 140-cm Frontrunner XP CTO catheter.
	Outback LTD Re-Entry Catheter	5.9	.014	120	.079	Highly visible "L" and "T" markers to easily reorient the re-entry cannula toward the true lumen.
Kensey Nash Corporation	Safe-Cross RF Total Occlusion Crossing Wires	n/a	.014, .035	175, 275	n/a	RF Energy for power and Optical Coherence Reflectometry for guidance in a wire-based platform.
Medtronic	Pioneer IVUS-Guided Re-Entry Catheter	6.9	.014	120	.092	64-phased array IVUS for accurate needle placement.

RADIOPAQUE MARKING TAPE

Company Name	Product Name	Working Length (cm)	Comments
LeMaitre Vascular	Glow 'N Tell Tape	30	Centimeter markings on flexible radiopaque ruler.
LeMaitre Vascular	LeMaitre Stent Guide	27	Millimeter markings on flexible radiopaque ruler.

VASCULAR PLUGS

Company Name	Product Name	Material Used	Sizes (mm)	Delivery Cable Length (cm)	Recommended Guide Catheter Sizes ID (inch)	Indicated Use
AGA Medical Corporation	Amplatzer Vascular Plug	Nitinol	4-8	135	.056	Arterial and venous embolizations in the peripheral vasculature.
			10-12	135	.067	
			14-16	135	.088	

TARGET RENAL THERAPY

Company Name	Product Name	Sheath Size (F)	System Lengths (cm)	Maximum Guiding Catheter Size (F)	Bifurcated Infusion Branch Profile (inch)	Indication
FlowMedica, Inc.	Benephit PV Infusion System	5	40, 105, 140	n/a	.035	Physician-specified agents to the peripheral vasculature including but not limited to the renal arteries.
	Benephit CV Infusion System	8	30, 35, 40, 45	6	.035	Physician-specified agents to the peripheral vasculature including but not limited to the renal arteries.