

OTHER DEVICES

SNARES					
Company Name	Product Name	Catheter Size/Length (F/cm)	Working Diameters (mm)	Snare Length (cm)	Snare Type
Angiotech	EN Snare Mini Snare System	3.2/150	2–4, 4–8	175	Three interlaced loops that self-size to the vessel diameter
	EN Snare Standard Snare System	6/100, 7/100	6–10, 9–15, 12–20, 18–30, 27–45	120	
Cook Medical	Curry Intravascular Retriever	8/100, 6.3/100	.021 inch, .018 inch	150	Folds at midpoint and forms a loop when passed through the catheter; supplied sterile in peel-open packages; intended for one-time use
	Dotter Intravascular Retriever	8/95	30-mm basket	7-cm basket	Helical loop basket is loaded into the hub of the catheter and expands upon emergence from the distal tip; supplied sterile in peel-open packages; intended for one-time use
ev3, Inc.	Amplatz Goose Neck Snare Kit	4/102	5, 10	120	Single 90° loop configuration
		16/162		65	
	Amplatz Goose Neck Snare Kit	6/102	15, 20, 25, 30, 35	120	
		16/223		65	
	Amplatz Goose Neck Microsnare Kit	2.3–3/150	2, 4, 7	175	
		2.3–3/175		200	
Radius Medical Technologies, Inc.	Expro Elite	.035 inch	5, 10, 15, 25, 35	150	Helical loop
	Micro Elite	.014 inch	2, 4, 7	180	

TARGETED RENAL THERAPY						
Company Name	Product Name	Sheath Size (F)	System Length (cm)	Maximum Guiding Catheter Size (F)	Bifurcated Infusion Branch Profile (inch)	Indicated Use
FlowMedica, Inc.	Benephit PV Infusion System	5	40, 105, 140	n/a	.035	Physician-specified agents to the peripheral vasculature including but not limited to the renal arteries
	Benephit XT Infusion System	5	35, 105			
	Benephit CV Infusion System	8	30, 35, 40, 45	6		