

AAA ENDOGRAFTS AND ACCESSORIES

MAIN BODY GRAFTS

| Company Name | Product Name | Main Body | | Iliac Leg(s) | |
|---------------------------|---------------------|--|----------------------------|----------------------------------|--|
| | | Length (cm) | Diameter (mm) | Length (cm) | Diameter (mm) |
| Cook Medical | Zenith Flex | 11.2, 12.5,* 12.6, 14.1, 14.3,* 15.5, 16.1,* 17, 17.9* | 22, 24, 26, 28, 30, 32, 36 | 3.7, 5.4, 7.1, 8.8, 10.5,† 12.2† | 8, 10, 12, 14, 16, 18, 20, 22, 24 |
| Endologix | Powerlink | 12, 13.5, 14, 15.5, 17.5+** | 25, 28, 34† | 4, 5.5 | 16 |
| Gore & Associates | Excluder | 12, 14, 15,† 16, 17,† 18 | 23, 26, 28.5, 31† | 9.5, 10, 11.5, 12, 13.5, 14 | 12, 14.5, 16, 18, 20 |
| Medtronic Cardio-Vascular | AneuRx AAAAdvantage | 13.5, 16.5 | 20, 22, 24, 26, 28 | 8.5, 11.5, 13.5 | 12, 13, 14, 15, 16, 18, 20 (20, 22, 24 distal flare) |

*Lengths for 36-mm-diameter grafts only

†Sizes available outside US

+Lengths for 8-, 10-, 12-mm-diameter grafts only

**Unibody main sizes 8, 10, 12

ANCILLARY COMPONENTS

| Company Name | Product Name | Main Body Extensions | |
|--------------------------|---------------------|--------------------------------|----------------------------|
| | | Length (cm) | Diameter (mm) |
| Cook Medical | Zenith | 3.9, 5, 5.8, 7.3 | 22, 24, 26, 28, 30, 32, 36 |
| | Renu | 4.3, 5.4,* 6.2, 7.7,* 8.1, 10* | 22, 24, 26, 28, 30, 32, 36 |
| Endologix | Powerlink | 5.5, 7.5, 8,† 9.5,† 10† | 25, 28, 34† |
| Gore & Associates | Excluder | 3.3, 4.5† | 23, 26, 28.5, 32† |
| Medtronic CardioVascular | AneuRx AAAAdvantage | 4 | 20, 22, 24, 26, 28 |

*Lengths for 36-mm-diameter grafts only

†Sizes available outside US

AAA ENDOGRAFTS AND ACCESSORIES

| Main Body Delivery System Profile | | | Fixation Location | Stent Expansion | Stent Material | Graft Material | Markets Available |
|-----------------------------------|-------------------------------|------------------------------------|---------------------------|-----------------|-----------------------|-----------------|--|
| Device Outer Diameter (F) | Sheath Required for Delivery? | Delivery Sheath Outer Diameter (F) | | | | | |
| 20, 23, 26 | No | n/a | Suprarenal | Self-expanding | Stainless steel | Woven polyester | US, Europe, Asia, Australia, Canada, South America, Middle East, Africa, Japan |
| 21 | No | n/a | Infrarenal or suprarenal† | Self-expanding | Cobalt-chromium alloy | ePTFE | US, Europe, Latin America, Africa, filed in Japan |
| 18, 20† | Yes | 20.5, 22.5† | Infrarenal | Self-expanding | Nitinol | ePTFE | US, Europe, Asia, Australia, Latin America |
| 21 | No | n/a | Infrarenal | Self-expanding | Nitinol | Woven polyester | US, Europe, Middle East, Africa, Latin America, Asia, Australia, Canada, India |

| Iliac Leg Extensions | | Aorto-Uni-Iliac/Converter | | Occluders |
|----------------------|--|--|----------------------------|----------------|
| Length (cm) | Diameter (mm) | Length (cm) | Proximal Diameter (mm) | Diameter (mm) |
| n/a | n/a | 8, 8.2 | 24, 28, 32, 36 | 14, 16, 20, 24 |
| n/a | n/a | 9.1, 9.9,* 10.8, 11.3, 11.6,* 12.5, 12.7,* 13.0, 13.3,* 14.4,* 14.7, 16.1* | 22, 24, 26, 28, 30, 32, 36 | n/a |
| 5.5, 8.8 | 16, 20 | n/a | n/a | n/a |
| 7, 9.5, 11.5, 13.5 | 10, 12, 14.5, 16, 18, 20 | n/a | n/a | n/a |
| 5.5, 8.5 | 12, 13, 14, 15, 16, 18, 20 (20, 22, 24 distal flare) | n/a | n/a | n/a |