

# BALLOONS

## ENDOVASCULAR TREATMENT AND MEASUREMENT

Company Name	Product Name	Measurement Range (mm)	Length (mm)	PTCA Balloon Diameter (mm)	PTCA Balloon Length (mm)	PTCA Nominal Pressure (atm)	
Medical Ventures Corp. (Angiometrx Inc.)	Metricath Gemini	2-4	6	2.5, 3, 3.5, 4	15	6	

Rated Burst Pressure (atm)	Catheter Usable Length (cm)	Maximum Guidewire (inch)	Minimum Guide Catheter Inner Diameter (inch)	Comments	FDA Indicated Use
16	140	.014 (.356 mm)	.064 (6 F)	Treatment balloon combined with a low-pressure measurement balloon that quickly and easily determines precise vessel and in-stent diameter and cross-sectional area	510(k) cleared for use in the renal, iliac, and femoral arteries