

# OTHER DEVICES

SNARES					
Company Name	Product Name	Catheter Size/ Length (F/cm)	Working Diameters (mm)	Snare Length (cm)	Snare Type
Abbott Vascular/ Vascular Solutions, Inc.	Expro Elite	.035 inch	5, 10, 15, 25, 35	150	Helical loop
Angiotech Interventional	EN Snare Mini Snare System	3.2/150	2-4, 4-8	175	Three interlaced nitinol loops that self-size to the vessel diameter
	Standard Snare System	6/100, 7/100	6-10, 9-15, 12-20, 18-30, 27-45	120	
Cook Medical	Curry Intravascular Retriever	8/100, 6.3/100	.021 inch, .018 inch	150	Folds at midpoint and forms a loop when passed through the catheter; supplied sterile in peel-open packages intended for one-time use
	Dotter Intravascular Retriever	8/95	30-mm basket	7-cm basket	Helical loop basket is loaded into the hub of the catheter and expands upon emergence from the distal tip; supplied sterile in peel-open packages intended for one-time use
ev3 Inc.	Amplatz Goose Neck Snare Kit	4/102	5, 10	120	Single 90° loop configuration
		16/162		65	
		6/102	15, 20, 25, 30, 35	120	
		16/223		65	
	Amplatz Goose Neck Microsnare Kit	2.3-3/150	2, 4, 7	175	
2.3-3/175		200			
Vascular Solutions, Inc.	Micro Elite Snare	.014 inch	2, 4, 7	180	Helical loop
	Expro Elite Snare	.035 inch	5, 10, 15, 25, 35	150	