

OCCCLUSION/MOLDING BALLOONS

| Company Name | Product Name | Size (F) | Recommended Introducer Size (F) | Wire Guide Diameter (inch) | Catheter Length (cm) | Maximum Inflated Balloon Diameter (mm) |
|-------------------------------|--------------------------------------|----------|---------------------------------|----------------------------|----------------------|--|
| Boston Scientific Corporation | Equalizer Balloon Catheter | 7 | 14 | .038 | 65, 100 | 20, 27, 33, 40 |
| Cook Incorporated | Coda Balloon Catheter | 10 | 14 | .035 | 100 | 32 |
| Medtronic | Reliant Stent Graft Balloon Catheter | 8 | 12 | .038 or smaller | 100 | 46 |