

TAA GRAFTS AND ACCESSORIES

TAA GRAFTS

Company Name	Product Name	Stent Material	Graft Material	Stent Expansion	Length (cm)	Diameter (mm)	Sheath Required for Delivery	Delivery Sheath (F)
Gore & Associates	TAG	Nitinol	ePTFE	Self-expanding	10, 15, 20	26, 28, 31, 34, 37, 40	Yes	20, 22, 24

WIRE GUIDES (BOSTON SCIENTIFIC CORPORATION)

Company Name	Product Name	Diameter (inch)	Length (cm)	Flexible Tip Length (cm)
Boston Scientific Corporation	Amplatz Super Stiff ST	.035, .038	145	3.5
			260	1
	Amplatz Super Stiff J	.035	145, 180, 260	3
		.038	145, 180	
	Amplatz Super Stiff	.035	145, 180, 260	6
		.038		
Meier Wire C Curve	.035	260, 300	10	
Meier Wire J Curve	.035	185	15	

WIRE GUIDES (COOK MEDICAL)

Company Name	Product Name	Diameter (inch)	Length (cm)	Primary Curve Radius (mm)	Secondary Curve Radius (mm)	Flexible Tip Length (cm)
Cook Medical	Lunderquist DC	.035	260, 300	75	15	4
	Lunderquist	.035	260, 300	0	0	4

OCCCLUSION/MOLDING BALLOONS

Company Name	Product Name	Size (F)	Recommended Introducer Size (F)	Wire Guide Diameter (inch)	Catheter Length (cm)	Maximum Inflated Balloon Diameter (mm)
Boston Scientific Corporation	Equalizer Balloon Catheter	7	14	.038	65, 100	20, 27, 33, 40
Cook Medical	Coda Balloon Catheter	10	14	.035	120	40
Medtronic CardioVascular	Reliant Stent Graft Balloon Catheter	8	12	.038 or smaller	100	46

DIAGNOSTIC CATHETERS

Company Name	Product Name	Type (Angiographic, Sizing, or Flush)	Size (F)	Recommended Introducer Size (F)	Wire Guide Diameter (inch)	Catheter Length (cm)
AngioDynamics, Inc.	Accu-Vu Sizing Catheter	Sizing	4, 5	4, 5	.035, .038	70, 90, 100
Cook Medical	Beacon Tip	Angiographic Kumpe	5	5	.035	100
	Aurous	Sizing pigtail	5	5	.035	100

DILATORS

Company Name	Product Name	Size (F)	Dilator Length (cm)	Wire Guide Compatibility (inch)
Cook Medical	Hydrophilic Dilator Set	20, 24	30	.035

PRESSURE SENSING

Company Name	Product Name	Type	Delivery Details	Implant Dimensions	MRI Compatibility	Stent Graft Compatibility	Patient Limitations	Comments
CardioMEMS, Inc.	EndoSure Wireless Pressure Sensor	Implantable pressure sensor	Working length: 35 cm; overall length: 67 cm; tether length: 110 cm; 17-F outer diameter; compatible with .038-inch guidewire	30 mm X 1.5 mm X 5 mm	Yes, per IFU	All commercially available	None, other than stent graft limitations	Wirelessly relays real-time pressure information from inside TAA sac